SPEED CROWN, for CEREC

ENGLISH INSTRUCTIONS FOR USE

I. Introduction

KATANA Zirconia ONE SPEED CROWN is a pre-sintered zirconia block, which is designed for CEREC milling systems. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia ONE SPEED CROWN has 14 shade variations: (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, NW). KATANA Zirconia ONE SPEED CROWN consists of 4 gradated shade layers.

The general clinical benefit of this product is to restore tooth function for the following INTENDED USE/ INDICATIONS FOR USE.

II. Intended use/ Indications for Use

KATANA Zirconia ONE SPEED CROWN is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, inlays, onlays and veneers).

III. Sintering Program

If "CEREC SpeedFire" is used as the sintering furnace, follow the Operating Instructions of the "CEREC SpeedFire" when sintering the block. If a sintering furnace other than the "CEREC SpeedFire" is used, review the sintering schedules shown below:

Sintering program 1

Temperature	Programming Rate	Holding Time
Room Temp. — 1550°C (2822°F)	10°C/min. (18°F/min.)	—
1550°C (2822°F)	—	2 hrs
1550°C (2822°F) — Room Temp.	-10°C/min. (-18°F/min.)	—

Sintering program 2

Temperature	Programming Rate	Holding Time
Room Temp. — 1560°C (2840°F)	35°C/min. (63°F/min.)	—
1560°C (2840°F)	—	0.5 hrs
1560°C (2840°F) — Room Temp. (*)	-45°C/min. (-81°F/min.)	—

*The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

Sintering program 3

Temperature	Programming Rate	Holding Time
Room Temp. — 900°C (1652°F)	130–342°C/min. (234–616°F/min.)	—
900°C (1652°F) — 1300°C (2372°F)	50–342°C/min. (90–616°F/min.)	—
1300°C (2372°F) — 1450°C (2642°F)	50–180°C/min. (90–324°F/min.)	—
1450°C (2642°F) — 1500°C (2732°F)	24°C/min. (43°F/min.)	—
1500°C (2732°F) — 1560°C (2840°F)	15–24°C/min. (27–43°F/min.)	—
1560°C (2840°F)	—	7–16 min.
1560°C (2840°F) — 1160°C (2120°F)	(-228)–(-70)°C/min. ((-410)–(-126)°F/min.)	—
1160°C (2120°F) — Room Temp. (*)	(-300)–(-70)°C/min. ((-540)–(-126)°F/min.)	—


*The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

Sintering program 4

Temperature	Programming Rate	Holding Time
Room Temp. — 1300°C (2372°F)	342°C/min. (616°F/min.)	—
1300°C (2372°F) — 1580°C (2876°F)	198°C/min. (356°F/min.)	—
1580°C (2876°F)	—	2 min.
1580°C (2876°F) — 1160°C (2120°F)	(-264)–(-70)°C/min. ((-475)–(-126)°F/min.)	—
1160°C (2120°F) — Room Temp. (*)	(-360)–(-70)°C/min. ((-648)–(-126)°F/min.)	—

*The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

IV. Image of product

This product has color gradation which varies from the enamel layer to the body layer, utilizing an additional two gradation layers between, for total of 4 layers in the a-axis direction.

The word "SPEED" printed on the block indicates that this product is "SPEED CROWN".

V. Composition

ZrO₂ (80-95%)
Y₂O₃ (5-15%)
Pigments (< 1%)

Units in parentheses are mass %.

VI. Type and Class (ISO 6872: 2015)

Type:II/ Class:4

VII. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9.8x10⁻⁶/K

VIII. Directions for Use

- Remove the block from the packaging and confirm that the block does not have a crack or other damage.
- Place the block into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restoration from the block and carefully smooth the surface where the restoration was connected to the mandrel with a diamond bur, etc.
- As desired, it can be polished the pre-sintered restoration with a fine diamond smoothing instrument, preferably silicone-free polishers.
- The cutting waste or dust, which is attached to the restoration, can be removed with a gentle air stream or a soft brush.
- Put the restoration into the refractory sagger tray and place them in the sintering furnace.
- All of the above processes are recommended to work under dry condition. In case of performing the process under wet conditions and using sintering furnace other than SpeedFire, dry the block at 200°C, 10 minutes in sintering furnace.
- Depending on the performance of the sintering furnace used, review the sintering schedule shown above (III. Sintering Program) before sintering the restorations.
- After sintering, adjust the restorations with a diamond bur as needed.
- Confirm that the restoration has no cracks.
- 1- FCZ:
 - Glaze baking: Create a high shine surface by polishing, especially on the contact areas, then apply the glaze on all surfaces in the usual manner.
 - Hand polishing: Create a high shine surface on the entire restoration by polishing without using the glaze. When the restoration is to be finished without the glaze, it is recommended to select one shade brighter than the final shade.
- 2 Frameworks: Build-up the dental porcelain (CERABIEN MiLai, CERABIEN ZR or CZR PRESS LF etc.) on the frameworks following the manufacturer's technical instructions. Check the coefficient of thermal expansion of the porcelain in the manufacturer's technical instructions to confirm compatibility.
- Use a cement that claims to be usable for zirconia. It is recommended to use resin cement to bond the restoration.

IX. Remarks on Handling

Contraindication:

- If the patient is hypersensitive to zirconia or any other components, this product must not be used.

Warning:

If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

Caution:

- This product should not be used when malocclusion, clenching or bruxism conditions are present.
- When milling the block or cutting, grinding and polishing the restoration, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.

- When milling the block or cutting, grinding and polishing the restoration, use safety glasses to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restorations. This product is only for dental application.
- The use of this product is restricted to dental professionals.
- Do not touch the items heated in the furnace with your bare hands.
- This product is particularly translucent; therefore, one should carefully consider its appropriateness for use when masking the underlying surface, such as a metal abutment, or a very dark or discolored tooth preparation.
- The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.
- When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.
- Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Veneer	0.4 mm or more*
Anterior crown or bridge	0.8 mm or more
Posterior crown, bridge	0.8 mm or more
Inlay and onlay	1.0 mm or more

* 0.4 mm or more of this product is for full zirconia veneers. Keep thickness 0.8 mm or more, if it is used for combination with the porcelain.

- Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations.
- When sintering this product in SpeedFire for the first time, confirm with the supplier of CEREC SpeedFire that it is compatible with this product.
- Check the coefficient of thermal expansion of the porcelain glaze in the manufacturer's technical instructions to confirm compatibility.
- Do not take the restorations out of the sintering furnace during high temperature, as the quenching causes the breaking. However, if an automatic opening type furnace is used under the sintering program 2-4, the restorations may be removed from the furnace at 800°C (1472°F) or less. When removed from the furnace, the restorations must be put on a tray made of ceramic fiber (for example: Noritake Porcelain Mat) to cool slowly.
- There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2-4. Please use the furnace after checking with the manufacturer that the furnace and temperature listed in sintering program 2-4 is compatible. If SpeedFire is used as a sintering furnace, remove the restorations from SpeedFire according to the guidance of sintering program.
- When using a sintering furnace for the first time and changing a sintering condition, colors after sintering may vary. Sinter a small piece of Zirconia beforehand and confirm the color.
- Only use this product to produce restorations for anterior crowns, posterior crowns, inlays/ onlays and veneers.

Caution in conjunction with usage:

- The carved seal side of "inLab" is the body layer; the opposite side is the enamel layer. Depending on the setting of the scanner or milling machine, the enamel layer and the body layer may be machined in the opposite direction. Confirm that the block is compatible before processing.
 - Dry milling is recommended. If wet milling/grinding is performed by using cooling water contaminated by silica-based glass ceramics (lithium disilicate glass, etc.), the translucency of the zirconia may be reduced after baking. Before wet milling/grinding, clean the milling/grinding chamber, cooling water tank and filter insert. The cooling water must be changed in order to assure optimum results.
 - Do not use the block if there is a crack noticed after removing it from the package.
 - If cutting waste or dust remains on the surface of the prosthesis, it may affect the color tone after sintering. When smoothing the surface of the prosthesis before sintering, preferably use a silicone-free polisher and/or ensure that cutting waste or dust is removed.
 - When smoothing the surface of the prosthesis before sintering, do not thin walls beyond the minimum wall thickness.
 - If a polished restoration will appear slightly darker compared to a glazed restoration of the same shade. If necessary, it is recommended to select a block shade that is one shade brighter than the target shade and adjust using stain.
 - If you find a crack in the restoration, do not use.
 - Dispose of this product as a medical waste to prevent infection.
- Storage:**
- Store in a cool and dry place. Keep away from direct sunlight.
 - The product should be stored at 10-30°C (50-86°F).
 - Do not remove the block from its packaging during storage.
 - The block is fragile, and requires care when handling.
 - The product must be stored in an appropriate place where only dental personnel have access.
 - The product must be used by the expiration date indicated on the package.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

"KATANA" and "CERABIEN" are registered trademarks or trademarks of NORITAKE CO., LIMITED.

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US Distributed by **Kuraray America, Inc.**
Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a dental professional.