# For 12Z and 14Z

(D1-025)  
Dental Zirconia

## For CEREC STML (Super Translucent Multi Layered) & ST (Super Translucent)

### ENGLISH INSTRUCTIONS FOR USE

**I. Introduction**  
KATANA Zirconia Block is a pre-sintered zirconia block which is designed for CEREC milling systems. (Please refer to your milling system's technical instructions regarding KATANA Zirconia Block.) KATANA Zirconia Block STML has 2 sizes: 12Z and 14Z and has 14 shade variations (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, NW). KATANA Zirconia Block ST has a 1 size: 12Z and 14Z and has 11 shades (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, NW). KATANA Zirconia Block STML consists of 4 graded shade layers; whereas KATANA Zirconia Block ST is mono color.

**II. Intended Use**  
KATANA Zirconia Block is used for the fabrication of the ceramic restorations (FCZ crowns, inlays, onlays and veneers).

**III. Sintering Program**  
If "CEREC SpeedFire" is used as the sintering furnace, follow the Operating Instructions of the "CEREC SpeedFire" when sintering the block. If a sintering furnace other than the "CEREC SpeedFire" is used, review the sintering schedules shown below.

Temperature	Programming Rate	Holding Time
Room Temp. — 1550°C (2822°F)	10°C/min. (18°F/min.)	—
1550°C (2822°F) — Room Temp.	-10°C/min. (-18°F/min.)	2 hrs

\*The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

Temperature	Programming Rate	Holding Time
Room Temp. — 1560°C (2840°F)	35°C/min. (63°F/min.)	—
1560°C (2840°F) — Room Temp. (+)	-45°C/min. (-81°F/min.)	0.5 hrs
1560°C (2840°F) — Room Temp. (+)	-45°C/min. (-81°F/min.)	—

\*The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

Temperature Programming Rate Holding Time  
Room Temp. — 900°C (1652°F) 130°C/min. (234°F/min.) —  
900°C (1652°F) — 1500°C (2732°F) 50°C/min. (90°F/min.) —  
1500°C (2732°F) — 1560°C (2840°F) 15°C/min. (27°F/min.) —  
1560°C (2840°F) — Room Temp. (+) -70°C/min. (-126°F/min.) —  
16 min  
The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

**IV. Image of product**  
Block Block (Body layer) Metallic holder  
112z x 15.3 x 19.2 x ±0.20mm  
142z x 17.8 x 19.2 x ±0.20mm  
16 mm

The STML is each color gradation which varies from the enamel layer to the body layer, utilizing an additional two gradation layers between, for total of 4 layers in the a-axis direction.

**V. Composition**  
ZrO<sub>2</sub> / Y<sub>2</sub>O<sub>3</sub>  
**VI. Type and Class**(ISO6872:2015)  
Type II / Class 4  
**VII. Physical Properties**  
Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9.8x10<sup>-6</sup>/°C  
**VIII. Directions for Use**  
(1) Remove the block from the packaging and confirm that the block does not have a crack or other damage.  
(2) Place the block into the milling machine, then begin the milling process following the milling system's technical instructions.  
(3) After milling, remove the restoration from the block with a diamond bur, etc.  
(4) The cutting waste or dust, which is attached to the restoration, can be removed with a gentle air stream or a soft brush.

Put the restoration into the refractory sagger tray and place them in the sintering furnace.  
(9) After sintering, adjust the restorations with water under dry condition.  
(10) After sintering, adjust the restorations with a diamond bur as needed.  
(11) Hand polishing: Create a high shine on the surface of the restoration by polishing without using a glaze. When the restoration is to be finished without the glaze, it is recommended to select one shade brighter than the intended color.  
(12) Use a resin cement to bond the restoration.  
**II. Remarks on Handling**  
1. If the patient is hypersensitive to zirconia or any other components, this product must not be used.  
2. CEREC SpeedFire is not recommended for anterior crowns, posterior crowns, inlays/onlays and veneers.  
**Warning:**  
(1) CEREC SpeedFire is a high speed sintering furnace. If used as intended, it will produce restorations for anterior crowns, posterior crowns, inlays/onlays and veneers.  
**Caution:**  
1. This product should not be used when malocclusion, clenching or bruxism conditions are present.  
2. Remove the block after cutting, grinding and polishing the restoration, use an approved dust mask and vacuuming with a filter to protect your lungs from inhaled dust.  
3. When milling the block and cutting, grinding and polishing the restoration, use safety glasses to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.  
4. Do not use for any purposes except for dental restorations. This product is for dental application.  
5. CEREC SpeedFire is designed to be used with the furnace with your brand. This product is particularly translucent; therefore, one should carefully consider its appropriateness for use when making the underlying surface, such as a metal abutment, or a very dark or discolored tooth preparation.  
6. The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.  
7. When preparing teeth, avoid the following: deep shoulders; J-margins; guide grooves; serrated margins; non-lapped abutments; undercut; glaze, glue, grooves; and sharp corners.  
8. Keep the following thickness of this product for fabricating prosthetics:  

Location & Indication	Wall thickness
Anterior crown	0.8 mm or more
Anterior crown	1.0 mm or more
Posterior crown, inlay and onlay	1.0 mm or more

  
9. Choose a shade color that is brighter than the intended color for a thick restoration. It may be used with the restoration furnace.  
10. When sintering this product in SpeedFire for the first time, confirm with the supplier of CEREC SpeedFire that it is compatible with this product.  
11. Check the coefficient of thermal expansion in the packaging as stated in the manufacturer's technical instructions to confirm compatibility.  
12. Do not take the restorations out of the sintering furnace during high temperature. If the restorations are removed during cooling, they may be deformed.  
13. An automatic opening type furnace is used under the sintering program 2 or 3, the restorations may be removed from the furnace at 800°C (1472°F) or less. When the restorations are removed from the furnace, they may be cooled in a tray made of ceramic fiber (for example, NortonKor Porcelain Mat) to cool slowly.  
14. There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2 or 3. Please use the furnace after consulting with the manufacturer that the furnace and temperature listed in sintering program 2 or 3 is compatible. If SpeedFire is used as a sintering furnace, remove the restorations from SpeedFire according to the guidance of sintering program.  
**(WARNING)**  
Kuraray NortonKor Dental Inc. will replace any product that is proven to be defective. Kuraray NortonKor Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.  
**(NOTE)**  
"KATANA" is trademarks of NORITAKE CO., LIMITED.

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